

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CHAC1 Polyclonal Antibody

多克隆抗体

产品货号: PA09068

第三版

**描述:**CHAC1 (ChaC, cation transport regulator homolog 1 (E. coli)) is a 264 amino acid cytoplasmic protein that belongs to the ChaC family. Induced by chemical activators of the unfolded protein response (UPR) pathway, including DTT, tunicamycin and thapsigargin, CHAC1 is a pro-apoptotic member of this pathway and may mediate the pro-apoptotic effects of the ATF4-ATF3-DDIT3/CHOP cascade. Overexpression of CHAC1 may lead to an increase in apoptosis, while suppression of CHAC1 results in subdued apoptosis.

**抗原:**Recombinant protein of human CHAC1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**