

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human CIAPIN1 Polyclonal Antibody

多克隆抗体

产品货号: PA09078

---

第三版

**描述:**Anamorsin is a protein that in humans is encoded by the CIAPIN1 gene. CIAPIN1 is a cytokine-induced inhibitor of apoptosis with no relation to apoptosis regulatory molecules of the BCL2 or CASP families. Expression of CIAPIN1 is dependent on growth factor stimulation.

**抗原:**Recombinant protein of human CIAPIN1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Ubiquitously expressed. Highly expressed in heart, liver and pancreas.