

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CKAP2 Polyclonal Antibody

多克隆抗体

产品货号: PA09083

第三版

**描述:**This gene encodes a cytoskeleton-associated protein that stabilizes microtubules and plays a role in the regulation of cell division. The encoded protein is itself regulated through phosphorylation at multiple serine and threonine residues. There is a pseudogene of this gene on chromosome 14. Alternative splicing results in multiple transcript variations.

**抗原:**Recombinant protein of human CKAP2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Abundant in testis, thymus, and in tumor derived cell lines, while barely detectable in liver, prostate, and kidney.