

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human HLTF Polyclonal Antibody

多克隆抗体

产品货号: PA09096

第三版

描述: This gene encodes a member of the SWI/SNF family. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein contains a RING finger DNA binding motif. Two transcript variants encoding the same protein have been found for this gene. However, use of an alternative translation start site produces an isoform that is truncated at the N-terminus compared to the full-length protein.

抗原: Synthetic peptide of human HLTF

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in brain, heart, kidney, liver, lung, pancreas, placenta and skeletal muscle.