

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GPR172B Polyclonal Antibody

多克隆抗体

产品货号: PA09125

第三版

**描述:**Biological redox reactions require electron donors and acceptor. Vitamin B2 is the source for the flavin in flavin adenine dinucleotide (FAD) and flavin mononucleotide (FMN) which are common redox reagents. This gene encodes a member of the riboflavin (vitamin B2) transporter family. Haploinsufficiency of this protein can cause maternal riboflavin deficiency. Multiple alternatively spliced variants, encoding the same protein, have been identified.

**抗原:**Synthetic peptide of human SLC52A1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Widely expressed. Highly expressed in the testis, placenta and small intestine. Expressed at lower level in other tissues.