

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human c-Rel Polyclonal Antibody

多克隆抗体

产品货号: PA09127

第三版

**描述:**The REL gene encodes c-Rel, a transcription factor that is a member of the Rel/NF $\kappa$ B family, which also includes RELA (MIM 164014), RELB (604758), NF $\kappa$ B1 (MIM 164011), and NF $\kappa$ B2 (MIM 164012). These proteins are related through a highly conserved N-terminal region termed the 'Rel domain,' which is responsible for DNA binding, dimerization, nuclear localization, and binding to the NF $\kappa$ B inhibitor

**抗原:**Recombinant protein of human REL

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**