

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CUEDC2 Polyclonal Antibody

多克隆抗体

产品货号: PA09149

第三版

描述:CUEDC2 is involved in ubiquitin- and proteasome-mediated degradation of progesterone receptor and estrogen receptor (ER)-alpha. CUEDC2 reduced PRB protein content and promoted progesterone-induced PRB degradation via the ubiquitin-proteasome pathway.

抗原:Recombinant protein of human CUEDC2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Significantly up-regulated in breast tumor tissues compared with matched adjacent normal tissues (at protein level). Levels inversely correlate with ESR1 in breast cancers and are lower in low-grade tumors compared to high-grade tumors.