

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GLMN Polyclonal Antibody

多克隆抗体

产品货号: PA09161

第三版

**描述:** This gene encodes a phosphorylated protein that is a member of a Skp1-Cullin-F-box-like complex. The protein is essential for normal development of the vasculature and mutations in this gene have been associated with glomuvenous malformations, also called glomangiomas. Alternatively spliced variants that encode different protein isoforms have been described but the full-length nature of only one has been determined.

**抗原:** Synthetic peptide of human GLMN

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Ubiquitous.