

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APOBEC3C Polyclonal Antibody

多克隆抗体

产品货号: PA09191

第三版

描述: This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control.

抗原: Synthetic peptide of human APOBEC3C

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expressed in spleen, testes, peripheral blood lymphocytes, heart, thymus, prostate and ovary.