

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SHROOM2 Polyclonal Antibody

多克隆抗体

产品货号: PA09197

第三版

**描述:**The protein encoded by this gene shares significant similarities with the apical protein from *Xenopus laevis* which is implicated in amiloride-sensitive sodium channel activity. This gene is a strong candidate gene for ocular albinism type 1 syndrome.

**抗原:**Synthetic peptide of human SHROOM2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Abundant in retina and melanoma; also in brain, placenta, lung, kidney and pancreas.