

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FOLH1B Polyclonal Antibody

多克隆抗体

产品货号: PA09218

第三版

描述: PSM is an integral membrane protein, and PCTA-1 is related to the galectin gene family, which mediate both cell-cell and cell-matrix interactions in a manner similar to the selectin subgroup of C-type lectins. STEAP is a serpentine transmembrane cell-surface tumor-antigen that is predicted to functions as a channel or transporter protein. In addition to prostate cancers, STEAP is also upregulated in bladder, colon and ovarian cancers.

抗原: Synthetic peptide of human FOLH1B

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Kidney and liver. Not expressed in the prostate.