

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PIWIL3 Polyclonal Antibody

多克隆抗体

产品货号: PA09255

第三版

描述: This gene encodes a member of the PIWI subfamily of Argonaute family proteins. This subfamily of proteins contains a PAZ domain, found in proteins involved in RNA-mediated gene silencing, and a C-terminal Piwi domain. The encoded protein is thought to function in maintenance of germline cells. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原: Synthetic peptide of human PIWIL3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed in testis.