

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse AMER1 Polyclonal Antibody

多克隆抗体

产品货号: PA09269

第三版

**描述:**The protein encoded by this gene upregulates transcriptional activation by the Wilms tumor protein and interacts with Human, Mouse and other proteins, including CTNNB1, APC, AXIN1, and AXIN2. Defects in this gene are a cause of osteopathia striata with cranial sclerosis (OSCS).

**抗原:**Synthetic peptide of human AMER1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Detected in fetal and adult kidney, brain and spleen.1 Publication