
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse DHRS3 Polyclonal Antibody

多克隆抗体

产品货号: PA09289

第三版

描述: Short-chain dehydrogenases/reductases (SDRs), such as DHRS3, catalyze the oxidation/reduction of a wide range of substrates, including retinoids and steroids.

抗原: Recombinant protein of human DHRS3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed with highest levels found in heart, placenta, lung, liver, kidney, pancreas, thyroid, testis, stomach, trachea and spinal cord. Lower levels found in skeletal muscle, intestine and lymph node. No expression detected in brain. In the retina, expressed in cone but not rod outer segments.