

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse DUSP23 Polyclonal Antibody

多克隆抗体

产品货号: PA09320

---

第三版

**描述:**Dual specificity protein phosphatase 23, also known as low molecular mass dual specificity phosphatase 3 (LDP-3), is an enzyme that in humans is encoded by the DUSP23 gene.

**抗原:**Recombinant protein of human DUSP23

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Widely expressed. Highly expressed in spleen, prostate, colon, adrenal gland, mammary gland, thyroid and trachea. Expressed at lower level in uterus, small intestine, bladder, bone marrow, brain, spinal cord and stomach.