

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ENDOG Polyclonal Antibody

多克隆抗体

产品货号: PA09321

第三版

描述:The protein encoded by this gene is a nuclear encoded endonuclease that is localized in the mitochondrion. The encoded protein is widely distributed among animals and cleaves DNA at GC tracts. This protein is capable of generating the RNA primers required by DNA polymerase gamma to initiate replication of mitochondrial DNA.

抗原:Synthetic peptide of human ENDOG

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: