

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CALCA Polyclonal Antibody

多克隆抗体

产品货号: PA09322

第三版

描述: This gene encodes the peptide hormones calcitonin, calcitonin gene-related peptide and katecalcitonin by tissue-specific alternative RNA splicing of the gene transcripts and cleavage of inactive precursor proteins. Calcitonin is involved in calcium regulation and acts to regulate phosphorus metabolism. Calcitonin gene-related peptide functions as a vasodilator while katecalcitonin is a calcium-lowering peptide. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原: Synthetic peptide of human CALCA

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: