

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human STK17A Polyclonal Antibody

多克隆抗体

产品货号: PA09351

第三版

描述: This gene is a member of the DAP kinase-related apoptosis-inducing protein kinase family and encodes an autophosphorylated nuclear protein with a protein kinase domain. The protein has apoptosis-inducing activity.

抗原: Synthetic peptide of human STK17A

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in placenta. Lower levels in heart, lung, skeletal muscle, kidney and pancreas.