

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EDN2 Polyclonal Antibody

多克隆抗体

产品货号: PA09405

第三版

**描述:** This gene encodes a member of the endothelin protein family of secretory vasoconstrictive peptides. The preproprotein is processed to a short mature form which functions as a ligand for the endothelin receptors that initiate intracellular signaling events. This gene product is involved in a wide range of biological processes, such as hypertension and ovulation.

**抗原:** Recombinant protein of human EDN2

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Expressed in lung, but not in placental stem villi vessels or cultured placental villi smooth muscle cells.