

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse PLA2G4B Polyclonal Antibody

多克隆抗体

产品货号: PA09464

第三版

描述: This gene encodes a member of the cytosolic phospholipase A2 protein family. Phospholipase A2 enzymes hydrolyze the sn-2 bond of phospholipids, releasing lysophospholipids and fatty acids. This enzyme may be associated with Human, Mouse mitochondria and early endosomes. Most tissues also express read-through transcripts from the upstream gene into this gene, some of which Human, Mouse encode fusion proteins combining the N-terminus of the upstream gene including its JmjC domain with the almost complete coding region of this gene, including the C2 and cytoplasmic phospholipase A2 domains.

抗原: Synthetic peptide of human PLA2G4B

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Widely expressed. Expressed at higher level in brain, heart, liver, cerebellum and pancreas. Isoform 3 is widely expressed.