



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FETUB Polyclonal Antibody

多克隆抗体

产品货号: PA09491

第三版

描述:The protein encoded by this gene is a member of the fetuin family, part of the cystatin superfamily of cysteine protease inhibitors. Fetuins have been implicated in several diverse functions, including osteogenesis and bone resorption, regulation of the insulin and hepatocyte growth factor receptors, and response to systemic inflammation. This protein may be secreted by cells.

抗原:Recombinant protein of human FETUB

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Liver and testis.