

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TSPYL2 Polyclonal Antibody

多克隆抗体

产品货号: PA09495

第三版

描述: This gene encodes a member of the testis-specific protein Y-encoded, TSPY-like/SET/nucleosome assembly protein-1 superfamily. The encoded protein is localized to the nucleolus where it functions in chromatin remodeling and as an inhibitor of cell-cycle progression. This protein may play a role in the suppression of tumor growth.

抗原: Synthetic peptide of human TSPYL2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitously expressed, with highest levels in brain, testis and heart, and lowest levels in liver and pancreas.