

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse GABPB2 Polyclonal Antibody

多克隆抗体

产品货号: PA09556

---

第三版

**描述:**The GABPB gene encodes at least 4 beta subunit isoforms that were originally thought to be encoded by 2 separate genes, GABPB1 and GABPB2. The gene on Hsa1q was named GABPB2 and the gene on Hsa15q was named GABPB1. The predicted GABPB1 and GABPB2 subunits contain 383 and 347 amino acids, respectively, and are highly homologous.

**抗原:**Recombinant protein of human GABPB2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**