

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NOD2 Polyclonal Antibody

多克隆抗体

产品货号: PA09557

第三版

描述: This gene is a member of the Nod1/Apaf-1 family and encodes a protein with two caspase recruitment (CARD) domains and six leucine-rich repeats (LRRs). The protein is primarily expressed in the peripheral blood leukocytes. It plays a role in the immune response to intracellular bacterial lipopolysaccharides (LPS) by recognizing the muramyl dipeptide (MDP) derived from them and activating the NF κ B protein. Mutations in this gene have been associated with Crohn disease and Blau syndrome.

抗原: Synthetic peptide of human NOD2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Monocytes-specific.