

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CARD17 Polyclonal Antibody

多克隆抗体

产品货号: PA09607

第三版

描述:Regulator of procaspase-1/CASP1 activation implicated in the regulation of the proteolytic maturation of pro-IL-1beta/IL1B and its release during inflammation. Inhibits the release of IL1B in response to LPS in monocytes. However, unlike CASP1, do not induce NF-kappa-B activation.

抗原:Synthetic peptide of human CARD17

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitous.