

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NCAPG2 Polyclonal Antibody

多克隆抗体

产品货号: PA09616

第三版

描述: This gene encodes a protein that belongs to the Condensin2nSMC family of proteins. The encoded protein is a regulatory subunit of the condensin II complex which, along with the condensin I complex, plays a role in chromosome assembly and segregation during mitosis. A similar protein in mouse is required for early development of the embryo. Alternate splicing results in multiple transcript variants.

抗原: Synthetic peptide of human NCAPG2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: