

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AGAP2 Polyclonal Antibody

多克隆抗体

产品货号: PA09665

第三版

**描述:**The protein encoded by this gene belongs to the centaurin gamma-like family. It mediates anti-apoptotic effects of nerve growth factor by activating nuclear phosphoinositide 3-kinase. It is overexpressed in cancer cells, and promotes cancer cell invasion. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

**抗原:**Synthetic peptide of human AGAP2

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Isoform 1 is brain-specific. Isoform 2 is ubiquitously expressed, with highest levels in brain and heart.