

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TSPYL2 Polyclonal Antibody

多克隆抗体

产品货号: PA09680

第三版

**描述:** This gene encodes a member of the testis-specific protein Y-encoded, TSPY-like/SET/nucleosome assembly protein-1 superfamily. The encoded protein is localized to the nucleolus where it functions in chromatin remodeling and as an inhibitor of cell-cycle progression. This protein may play a role in the suppression of tumor growth.

**抗原:** Synthetic peptide of human TSPYL2

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Ubiquitously expressed, with highest levels in brain, testis and heart, and lowest levels in liver and pancreas.