

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD80 Polyclonal Antibody

多克隆抗体

产品货号: PA09685

第三版

描述:The protein encoded by this gene is a membrane receptor that is activated by the binding of CD28 or CTLA-4. The activated protein induces T-cell proliferation and cytokine production. This protein can act as a receptor for adenovirus subgroup B and may play a role in lupus neuropathy.

抗原:Synthetic peptide of human CD80

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed on activated B-cells, macrophages and dendritic cells.