

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CRHR1 Polyclonal Antibody

多克隆抗体

产品货号: PA09718

第三版

**描述:** This gene encodes a G-protein coupled receptor that binds neuropeptides of the corticotropin releasing hormone family that are major regulators of the hypothalamic-pituitary-adrenal pathway. The encoded protein is essential for the activation of signal transduction pathways that regulate diverse physiological processes including stress, reproduction, immune response and obesity. Alternative splicing results in multiple transcript variants, one of which represents a read-through transcript with the neighboring gene MGC57346.

**抗原:** Synthetic peptide of human CRHR1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Predominantly expressed in the cerebellum, pituitary, cerebral cortex and olfactory lobe.