

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TDGF1 Polyclonal Antibody

多克隆抗体

产品货号: PA09727

第三版

描述: This gene encodes an epidermal growth factor-related protein that contains a cripto, FRL-1, and cryptic domain. The encoded protein is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. Mutations in this gene are associated with forebrain defects. Pseudogenes of this gene are found on chromosomes 2, 3, 6, 8, 19 and X. Alternate splicing results in multiple transcript variants.

抗原: Synthetic peptide of human TDGF1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Preferentially expressed in gastric and colorectal carcinomas than in their normal counterparts.