

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CALCRL Polyclonal Antibody

多克隆抗体

产品货号: PA09735

第三版

**描述:** Calcitonin receptor-like (CALCRL), also known as the calcitonin receptor-like receptor (CRLR), is a human protein. The protein encoded by the CALCRL gene is a G protein-coupled receptor related to the calcitonin receptor. CALCRL is linked to one of three single transmembrane domain receptor activity-modifying proteins (RAMPs) that are essential for functional activity

**抗原:** Synthetic peptide of human CALCRL

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$  g/ml.

**特异性:** Predominantly expressed in the lung and heart.