

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NPAP1 Polyclonal Antibody

多克隆抗体

产品货号: PA09822

第三版

**描述:** This gene is located in the Prader-Willi syndrome region on chromosome 15. This gene is biallelically expressed in adult testis and brain but is paternally imprinted in fetal brain. Defects in this gene may be associated with Prader-Willi syndrome.

**抗原:** Synthetic peptide of human NPAP1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Testis-specific in adults. In fetal brain expressed only from the paternal allele.