

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human KCNK17 Polyclonal Antibody

多克隆抗体

产品货号: PA09836

第三版

描述:The protein encoded by this gene belongs to the family of potassium channel proteins containing two pore-forming P domains. This channel is an open rectifier which primarily passes outward current under physiological K⁺ concentrations. This gene is activated at alkaline pH. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

抗原:Recombinant protein of human KCNK17

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: