

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EVC2 Polyclonal Antibody

多克隆抗体

产品货号: PA09892

第三版

描述: This gene encodes a protein that functions in bone formation and skeletal development. Mutations in this gene, as well as in a neighboring gene that lies in a head-to-head configuration, cause Ellis-van Creveld syndrome, an autosomal recessive skeletal dysplasia that is also known as chondroectodermal dysplasia. Mutations in this gene also cause acrofacial dysostosis Weyers type, also referred to as Curry-Hall syndrome, a disease that combines limb and facial abnormalities. Alternative splicing results in multiple transcript variants.

抗原: Synthetic peptide of human EVC2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Found in the heart, placenta, lung, liver, skeletal muscle, kidney and pancreas.1 Publication