

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ARHGAP4 Polyclonal Antibody

多克隆抗体

产品货号: PA09916

第三版

**描述:** This gene encodes a member of the rhoGAP family of proteins which play a role in the regulation of small GTP-binding proteins belonging to the RAS superfamily. The protein encoded by the orthologous gene in rat is localized to the Golgi complex and can redistribute to microtubules. The rat protein stimulates the activity of some Rho GTPases in vitro. Genomic deletions of this gene and a neighboring gene have been found in patients with nephrogenic diabetes insipidus. Multiple transcript variants encoding different isoforms have been found for this gene.

**抗原:** Synthetic peptide of human ARHGAP4

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Predominantly in hematopoietic cells (spleen, thymus and leukocytes); low levels in placenta, lung and various fetal tissues.