

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Chicken OVAL Polyclonal Antibody

多克隆抗体

产品货号: PA09955

第三版

描述: Ovalbumin (abbreviated OVA) is the main protein found in egg white, making up 60-65% of the total protein. Ovalbumin displays sequence and three-dimensional homology to the serpin superfamily, but unlike most serpins it is not a serine protease inhibitor. It is presumed to be a storage protein.

抗原: Recombinant protein of Chicken OVAL

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Major protein of egg white.