

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GAGE12I Polyclonal Antibody

多克隆抗体

产品货号: PA09966

第三版

描述: This gene belongs to a family of genes that are expressed in a variety of tumors but not in normal tissues, except for the testis. The sequences of the family members are highly related but differ by scattered nucleotide substitutions. The antigenic peptide YYWPRPRRY, which is also encoded by several other family members, is recognized by autologous cytolytic T lymphocytes.

抗原: Synthetic peptide of human GAGE7

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: