

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANGPT4 Polyclonal Antibody

多克隆抗体

产品货号: PA09971

第三版

描述:Angiopoietins are proteins with important roles in vascular development and angiogenesis. All angiopoietins bind with similar affinity to an endothelial cell-specific tyrosine-protein kinase receptor. The mechanism by which they contribute to angiogenesis is thought to involve regulation of endothelial cell interactions with supporting perivascular cells. The protein encoded by this gene functions as an agonist and is an angiopoietin.

抗原:Synthetic peptide of human ANGPT4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in the lung with much lower levels found in other tissues.