

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GLYR1 Polyclonal Antibody

多克隆抗体

产品货号: PA10003

第三版

**描述:**GLYR1 encodes a deduced 553-amino acid protein with a molecular mass of 60.5 kD. GLYR1 contains 2 nuclear localization signals, an N-terminal PWWP domain and an AT-hook, and a C-terminal NAD-binding domain. Western blot, RT-PCR, and immunohistochemistry showed that GLYR1 localizes to the nucleus in a variety of cell lines, with robust expression in mouse macrophage RAW274.7 cells.

**抗原:**Synthetic peptide of human GLYR1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**