

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL3RB Polyclonal Antibody

多克隆抗体

产品货号: PA10005

第三版

描述:The protein encoded by this gene is the common beta chain of the high affinity receptor for IL-3, IL-5 and CSF. Defects in this gene have been reported to be associated with protein alveolar proteinosis (PAP).

抗原:Synthetic peptide of human CSF2RB

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: