

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse CMC4 Polyclonal Antibody

多克隆抗体

产品货号: PA10028

---

第三版

**描述:** This gene was identified by involvement in some t(X;14) translocations associated with Human, Mouse and T-cell proliferations. This region has a complex gene structure, with a common promoter and 5' exon spliced to two different sets of 3' exons that encode two different proteins. This gene represents the downstream 8 kDa protein that localizes to mitochondria.

**抗原:** Recombinant protein of human CMC4

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Expressed in many tissues with a relatively high level in skeletal muscle.