

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CK-10 Polyclonal Antibody

多克隆抗体

产品货号: PA10043

第三版

描述:This gene encodes a member of the type I (acidic) cytokeratin family, which belongs to the superfamily of intermediate filament (IF) proteins. Keratins are heteropolymeric structural proteins which form the intermediate filament. These filaments, along with actin microfilaments and microtubules, compose the cytoskeleton of epithelial cells. Mutations in this gene are associated with epidermolytic hyperkeratosis. This gene is located within a cluster of keratin family members on chromosome 17q21.

抗原:Synthetic peptide of human KRT10

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Seen in all suprabasal cell layers including stratum corneum.