
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRPA1 Polyclonal Antibody

多克隆抗体

产品货号: PA10093

第三版

描述:The structure of the protein encoded by this gene is highly related to both the protein ankyrin and transmembrane proteins. The specific function of this protein has not yet been determined; however, studies indicate the function may involve a role in signal transduction and growth control.

抗原:Synthetic peptide of human TRPA1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed at very low level. Expressed at very low level in human fibroblasts and at a moderate level in liposarcoma cells.