

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse KIRREL2 Polyclonal Antibody

多克隆抗体

产品货号: PA10122

第三版

描述: NEPH3 (nephrin-like protein 3), also referred to as NLG1, FILTRIN or KIRREL2 (kin of IRRE like 2), is a 708 amino acid protein that is a member of the immunoglobulin superfamily of cell adhesion molecules. NEPH3 consists of five extracellular Ig-like repeats, a transmembrane domain, several glycosylation sites and a cytoplasmic domain that has a stretch of nine conserved residues. NEPH3 localizes to the interpodocyte filtration slit in the kidney and to the Langerhans islet beta cells in the pancreas. NEPH3 is thought to be involved in the maintenance of the glomerular filtration barrier in the kidney and in pancreas development.

抗原: Recombinant protein of human KIRREL2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in beta-cells of the pancreatic islets.