

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TNN Polyclonal Antibody

多克隆抗体

产品货号: PA10139

第三版

描述: This protein is involved in neurite outgrowth and cell migration in hippocampal explants. It has three EGF-like domains, one fibrinogen C-terminal domain and nine fibronectin type III domains. Tenascins are extracellular matrix proteins present during the development of organisms as well as in pathological conditions. Tenascin-W, the fourth and last member of the tenascin family remains the least well-characterized one.

抗原: Synthetic peptide of human TNN

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Not detected in normal adult mammary tissues or brain but expressed in most breast tumors and brain tumors, such as glioblastomas, astrocytomas and oligodendrogliomas, tested (PubMed:17909022, PubMed:19884327). In brain tumors, detected around the endothelial cell layer of the blood vessels (PubMed:19884327).