
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NT5C1A Polyclonal Antibody

多克隆抗体

产品货号: PA10143

第三版

描述:Cytosolic nucleotidases, such as NT5C1A, dephosphorylate nucleoside monophosphates. Dephosphorylates the 5' and 2'(3')-phosphates of deoxyribonucleotides and has a broad substrate specificity. Helps to regulate adenosine levels in heart during ischemia and hypoxia.

抗原:Recombinant protein of human NT5C1A

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in skeletal muscle. Detected at intermediate levels in heart, brain, kidney and pancreas.