

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IFNA2 Polyclonal Antibody

多克隆抗体

产品货号: PA10148

第三版

**描述:** This gene is a member of the alpha interferon gene cluster on chromosome 9. The encoded protein is a cytokine produced in response to viral infection. Use of the recombinant form of this protein has been shown to be effective in reducing the symptoms and duration of the common cold.

**抗原:** Synthetic peptide of human IFNA2

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**