

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human OCIAD2 Polyclonal Antibody

多克隆抗体

产品货号: PA10163

第三版

描述:OCIAD2 (OCIA domain containing 2) is a 154 amino acid protein that localizes to endosomes and contains one OCIA domain. OCIAD2 exists as multiple alternative spliced isoforms and is encoded by a gene which localizes to human chromosome 4. Representing approximately 6% of the human genome, chromosome 4 contains nearly 900 genes, one of which is the Huntingtin gene, which is found to encode an expanded glutamine tract in cases of Huntington's disease.

抗原:Recombinant protein of human OCIAD2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: