

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PAGE1 Polyclonal Antibody

多克隆抗体

产品货号: PA10189

第三版

描述: This gene belongs to a family of genes that are expressed in a variety of tumors but not in normal tissues, except for the testis. Unlike the other gene family members, this gene does not encode an antigenic peptide. Nothing is presently known about the function of this protein.

抗原: Recombinant protein of human PAGE1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: